


We are endeavouring to simplify the provision of minor adaptations.

To do this we have produced a small catalogue of standard rails.

The rails are numbered in the catalogue (this number has no technical implication).

We hope these cover most types of requests.


In future when referring work to any technician or Hull City building services, to fit rails or equipment on behalf of Hull Social Services, please use the number in the catalogue to describe your request e.g. catalogue reference no 1 rail next to toilet.



12" (300mm) PLASTIC  
FLUTED GRAB RAILS  
REFERENCE NO( 1 )

Plastic fluted rails are only to be used indoors.  
(We have found in the past that if fitted out doors  
they will eventually crack).


Hull Social Services  
01/01/08



18" (450mm) PLASTIC  
FLUTED GRAB RAILS

REFERENCE NO ( 2 )

Plastic fluted rails are only to be used indoors.  
(We have found in the past that if fitted out doors  
they will eventually crack).




24" (600mm) PLASTIC

FLUTED GRAB RAIL

REFERENCE NO ( 3 )

Plastic fluted rails are only to be used indoors.  
(We have found in the past that if fitted out doors  
they will eventually crack).



18" (445mm) MILD STEEL


FLAT END GRAB RAIL

CRANKED

REFERENCE NO ( 4 )

Hull Social Services  
01/01/08






22" (560mm) MILD STEEL

FLAT END GRAB RAIL

CRANKED

REFERENCE NO ( 5 )

Hull Social Services  
01/01/08




18" (445mm) MILD STEEL

FLAT END GRAB RAIL

STRAIGHT

REFERENCE NO ( 6 )

Hull Social Services  
01/01/08



22" (560mm) MILD STEEL


FLAT END GRAB RAIL

STRAIGHT

REFERENCE NO ( 7 )

Hull Social Services  
01/01/08






13" (320mm) PLASTIC

DERBY GRAB HANDLE

REFERENCE NO ( 8 )

Plastic Derby rails are only to be used indoors.  
(We have found in the past that if fitted out doors  
they will eventually crack).

Hull Social Services  
01/01/08



28" (707mm) METAL

NEWEL RAIL

LEFT HAND

REFERENCE NO ( 9 )




28" (707mm) METAL

NEWEL RAIL

RIGHT HAND

REFERENCE NO ( 10 )





FLOOR FIXED

TOILET FRAME

AND SEAT

REFERENCE NO ( 11 )

Hull Social Services  
01/01/08



HALF SURREY TOILET RAIL

REFERENCE NO ( 12 )

Hull Social Services  
01/01/08





FLOOR FIXED TOILET

SURROUND

REFERENCE NO ( 13 )

Hull Social Services  
01/01/08



30 " (769mm)  
DEVON MK11  
FOLDING SUPPORT  
RAIL  
REFERENCE NO 14

Hull Social Services  
01/01/08



30'' (760mm)  
DEVON RAIL WITH  
FOLDING LEG  
REFERENCE NO ( 15 )

Hull Social Services  
01/01/08